



Maxson®


Maxson®

Safety Personal
Protective Equipment

Maxson®

Autoclavable Goggles

Item No.: ST011



Specifications:

- Autoclavable at 250°F/121°C (40 cycles at 15 min.)
- Maxson Clear PC Lens Coating with Anti-Fog, Anti-Scratch & Anti-Splash (AF, AS, AS)
- Gamma Sterilizable (14 cycles at 10-15 kGy for labs or 14 cycles at 20-25 kGy for hospitals)
- Lightweight ~70-80g
- CE Marked + Medical Grade + PPE EN 166 + EU 2016/425 + EN 167 + EN 168 + EN71 as required.

Maxson Lenses are crystal clear,
provide safety in eyewear and has anti-fog coating





Autoclavable Goggles

Eye Shield

Snap N' Shield

Face Shield

Junior Face Shield

PED Face Shield

Autoclavable	√	-	Frame	-	-	-
Antimicrobial	-	-	Frame	-	-	-
Replaceable Lens	-	√	√	-	-	-
Adult	√	√	√	√		
Junior	√				√	
PED (Pediatrics)						√
Droplets Splash Guard	√	√	√	√	√	√
Breathability	√	√	√	√	√	√
Cleanroom Clearance (ISO 8 / 100,000 Class)	√	-	-	-	-	-
Medical Grade Clearance	√	√	√	√	√	√
PPE Clearance (EN 166, EU 2016/425)	P	-	P	NB.2575	P	P
Biocompatibility Clearance (ISO10993)	√	√	√	√	√	√
Latex Free	√	√	√	√	√	√
CE & FDA Clearance	√	√	√	√	√	√
ISO 13485 Factory Clearance	√	√	√	√	√	√

Alleva Medical's Maxson™ PPE offer a variety of designs and features with sizes from Adults, Junior to Pediatrics (0 to 12 months). The lenses and shields have anti-fog coating with durable structure meeting PPE as well as medical standards. The Snap N' Shield features antimicrobial properties that are 97-99% effective with Escherichia Coli and Staphylococcus Aureus as per the ISO and ASTM requirements.

The frame for the Snap N' Shield and the Autoclavable Goggles are autoclavable at 121°C for high temperature sterilization requirement of the hospitals and laboratories. Laser ID, branding and Made in Hong Kong products can be ordered as required. Please see our full line offering with accessories below.

Goggles & Eye Shield

ST011	MI364	SM002	MI149	MI149
				
Maxon Autoclavable Goggles	Maxon Splash Safety Goggles	Maxon Junior Safety Goggles	Eye Shield	Eye Shield Frame
(W) 173 x (H) 105 mm		(W) 155 x (H) 75 mm	(W) 240 x (H) 70 mm	

Face Shield

MI383-0001	MI383-0002	MI383-0003	MI163-0001	MI213	MI256
					
Maxon 1-Day Snap N' Shield	Maxon 7-Day Snap N' Shield	Maxon Snap N' Shield Lens	Maxon Face Shield	Maxon Junior Face Shield	Maxon PED Face Shield
Shield Size: (W) 320 x (H) 220 mm				(W) 320 x (H) 200 mm	(W) 230 x (H) 180 mm

Maxson®

ALLEVA MEDICAL LTD.

Suite M-Q, 12th Floor, Kings Wing Plaza 2, 1 On Kwan St. Shek Mun, Shatin, N.T., Hong Kong
Tel: +852 3166 7200 Fax: +852 2355 7663

ALLEVA MEDICAL USA LLC

Suite 205, 2460 General Armistead Avenue, Norristown, PA 19403
Tel: 1 267 428 2487 info@allevamedical.com



www.allevamedical.com

All products are manufactured under US FDA, ISO 13485, and Japan PMDA registered facility of Alleva Medical Ltd. Alleva Medical USA LLC, Suite 205, 2460 General Armistead Avenue, Norristown, PA 19403

Alleva Medical®, extriCARE®, Vide®, Maxson® are Trademarks and/or Registered Trademarks manufactured by Alleva Medical Ltd., its subsidiaries, and/or affiliates in the United States, Europe, China and other countries. Copyright 2021. All dimensions, weights, colours and contents may vary.