

to avoid mismatching





## **Enteral Feeding**

Item No.	Volume (ml)	Sterile	Packaging
MI191	1200	Non-Sterile	30pcs/carton
ST295	1200	EO	30pcs/carton
MI220	600	Non-Sterile	30pcs/carton
ST225	600	EO	30pcs/carton

## **Oral Syringe**

Item No.	Volume (cc)	Sterile	Packaging
MI194	1	-	-
MI196	3	-	100pcs/inner; 2400pcs/outer
MI197	5	-	100pcs/inner; 2400pcs/outer
MI198	10	-	100pcs/inner; 1800pcs/outer

- Barrel and plunger are assembled and ready for use with clear graduations
- Tip cap is provided for recap if multiple use is required
- Bright purple plunger alerts clinician that syringe if for oral use
- Special tip of the plunger cannot fit with hypodermic needle
- Material does not affect medications or solutions to be used
- Non-sterile bulk pack at 100pcs/poly bag

AS020	20	EO	50pcs/inner; 600pcs/outer

Maxson Diposable Catheter Tip Syringe (with ENFit TM connector)



10