



# ST382 Dry Seal Chest Drainage Unit Distributor Sell Sheet

# ST382 Dry Seal Chest Drainage Unit



Vide®'s Dry Seal is a 2000mL sterile disposable three-chamber system that offers the reliability of a traditional water seal drain and the convenience of a water-less operation. With a proprietary dry suction control technology, suction can be adjusted easily and precisely with one hand. Featuring needleless sampling port and relief valves, the dry seal drain is built entirely with patient safety and ergonomics in mind.

**Indications for use:** pneumothorax, hemothorax, mediastinal shift and malignant effusions.

## Features:



### Colored Liquids For Water Seal

Pre-filled syringe stored safely at the back of the drain and is ready for use



### Intuitive Sliding Suction Control panel

- no need for finding and filling sterile water!



### Suction and Negative Pressure Indicators

Allow easy visualization of the status of the unit



### 7-Scale Air Leak Monitor

for more precise air leak quantification



### Negative Pressure Relief Valve

release automatically during prolonged episode of negative intrapleural pressure



### Y-Shape Sampling Port

with self-sampling needless feature, allowing access to fresh blood samples



**Vide®'s Calibrated** Dry Suction Chest Drain can maintain consistent suction control at varying negative pressure or flow rate from the source



**Knock-over nozzle** to ensure function is retained even after unit is accidentally knocked over



### Low Positive Pressure Relief Valve

release automatically at 2cmH<sub>2</sub>O



**8-Discrete Suction Levels** with Suction Control Indicator  
5cmH<sub>2</sub>O to 40cmH<sub>2</sub>O

## Accessories:

**HC063-0001**

Chest Drain Basket  
with Clamp & Steel Hanger



**HC063-0002**

Chest Drain Basket with Clamp



**ST380-1031**

Small-bore Adaptor (22-26 Fr)



## Shipping / Packaging information:

Item No.	Type	Volume(mL)	Packaging
ST382	Dry Seal Chest Drainage Unit	2000mL	5pcs / Carton;



### ALLEVA MEDICAL LTD.

Suite M-Q, 12th Floor, Kings Wing Plaza 2, 1 On Kwan St. Shek Mun, Shatin, N.T., Hong Kong  
Tel: +852 3166 7200 Fax: +852 2355 7663

### ALLEVA MEDICAL USA LLC

Suite 205, 2460 General Armistead Avenue, Norristown, PA 19403  
Tel: 1 267 428 2487 info@allevamedical.com



[www.allevamedical.com](http://www.allevamedical.com)

ALL products are manufactured under US FDA, ISO 13485, and Japan PMDA registered facility of Alleva Medical Ltd.  
Alleva Medical USA LLC, Suite 205, 2460 General Armistead Avenue, Norristown, PA 19403

Alleva Medical®, extriCARE®, Vide®, Maxson® are Trademarks and/or Registered Trademarks manufactured by Alleva Medical Ltd., its subsidiaries, and/or affiliates in the United States, Europe, China and other countries. Copyright 2021.

All dimensions, weights, colours and contents may vary.