




# SURGICAL SUCTION

Distributor Sell Sheet

# The Suction System - Bag Liners

Our distinguished Vide® suction liners come in different sizes and bag options. Each liners type was designed to withstand a 1.5m drop and 50kg pressure test - as accidents do happen. Our newest design-patented bag liners are color-coded to provide instant identification of bag volumes. Each liner is equipped with anti-reflux valve at the patient port and a high efficiency filter at vacuum end to provide maximum safety. Pre-gelled liner option is available for those who require solidifiers but prefer immediate installation during usage.



- ① Ortho or Tandem Port for Multiple Canisters Linkage
- ② 90° 5-in-1 Elbow
- ③ Anti-Reflux Valve
- ④ Self Shut-off Ultra-high 99.9% Ultra-high Efficiency Filter\*
- ⑤ Solidifier
- ⑥ Colour Coding
- ⑦ Air tight port caps
- ⑧ Handle
- ⑨ Fold Flat Design
- ⑩ ON/OFF Switch

- Product Standard:
- ✓ FDA & CE Cleared
  - ✓ Conforms to ISO 10079-3:2015
  - ✓ Manufactured in ISO 13485 Facility
  - ✓ \*Filter Tested by Nelson Laboratories Inc.

## Technical Specifications & Parameters:

	Bag Liners				Bag Liners with Solidifiers				Outer Canisters				Canisters		
Item No.	MI285-0008	MI286-0008	MI287-0008	MI288-0008	MI285-0015	MI286-0015	MI287-0015	MI288-0015	MI129-0013	MI358-0013	MI301-0013	MI302-0013	MI055	MI056	MI202
Volume (ml)	1000	1500	2000	3000	1000	1500	2000	3000	1000	1500	2000	3000	800	1200	1500
Solidifier Pack	-				2x10g	3x10g	4x10g	6x10g	-				-		
Max. Vacuum	1000kPa (750 mmHg)				1000kPa (750 mmHg)				1000kPa (750 mmHg)				1000kPa (750 mmHg)		
Material	Lid: PE Blend Liner: Proprietary PE/PA Blend				Lid: PE Blend Liner: Proprietary PE/PA Blend				Canister: PC				Polycarbonate		
Packaging (Pcs/carton)	50	50	50	50	50	50	50	50	20	20	20	12	48	48	48
Volume (cbm)	0.072				0.075				0.106	0.112	0.128	0.106	0.133	-	0.166
Carton Dimension (cm)	59 x 49 x 25				59 x 49 x 26				70 x 36 x 42	69 x 36 x 45	69 x 64 x 29	75 x 33 x 43	57 x 44 x 45	-	64 x 52 x 50
G.W. (kg/carton)	5.3	5.6	5.6	5.9	5.9	6.8	6.8	9.2	8.6	8.7	10.6	7.7	-	-	-
Colour Code	Yellow 	Violet 	Aqua 	Green 	Yellow 	Violet 	Aqua 	Green 	Yellow 	Violet 	Aqua 	Green 	-		
Drop Test Height	1.5m				1.5m				-				-		

## System Accessories:

### Trolley & Stands

MI177-0007	MI177-0002	HC013-0002	HC025-0002	HC007
 1pc/Carton; 0.070CBM (4X Canister) 106cm	 1pc/Carton; 0.070CBM 106cm	 1pc/Carton; 0.038CBM 56cm	 1pc/Carton 37cm	 1pc/Carton 100cm
Regulator w/Metal Roller Stand	Metal Roller Stand	Metal Mini Stand	Tabletop Stand	Tabletop Stand








### Tubing & Connectors

HC005	ST258-0001	ST258-0002	HC033	ST001	MI168	S11-MI285-00020000
 50pc/Carton; 0.050CBM	 50pc/Carton; 0.050CBM	 50pc/Carton; 0.050CBM	 50pc/Carton; 0.050CBM	 50pc/Carton; 0.050CBM	 100pcs/Carton; 0.051CBM	
Roller Stand Hose	1.8m Connection Tube EO Sterile	3.6m Connection Tube EO Sterile	Tandem Tube	Suction Catheter EO Sterile	Connection Tube	Elbow Connector





## Wall Mounts & Brackets

MI167	HC051	HC034-0001	HC034-0002	MI129-0012	MI358-0012	MI301-0012	MI302-0012
200pcs/Carton; 0.032CBM 		1/1.5/2L 	3L 				
Hanger Plate	Hanger Adaptor (SR)	Hanger Adaptor (AS)	Hanger Adaptor (AS)	1L Canister	1.5L Canister	2L Canister	3L Canister

## Specialized Applications

MI153-00020000	MI153-00040000	MI153-00010000	HC026	MI291	MI359	MI361
750pcs/Inner; 4500pcs/Carton; 0.025CBM 		50pcs/Inner; 400pcs/Carton; 0.02CBM 				
2g Solidifier	10g Solidifier	20g Solidifier	Spicment Collector	15L Suction Canister	300ml Measurement Cup	250ml Suction Canister

## Yankauer Suction Tube

MI303* / ST303**	MI305* / ST305**	MI306* / ST306**	MI307* / ST307**
			
Square tip, suction control	Square tip, no control	Crown tip, suction control	Crown tip, no control

\*Non-sterile \*\*Sterile



### ALLEVA MEDICAL LTD.

Suite M-Q, 12th Floor, Kings Wing Plaza 2, 1 On Kwan St. Shek Mun, Shatin, N.T., Hong Kong  
Tel: +852 3166 7200 Fax: +852 2355 7663

### ALLEVA MEDICAL Devices

1050 W. Central Avenue, Unit B, Brea, CA 92821  
Tel: 1 267 428 2487 info@allevamedical.com



www.allevamedical.com

ALL products are manufactured under US FDA, ISO 13485, and Japan PMDA registered facility of Alleva Medical Ltd. Alleva Medical Devices, 1050 W. Central Avenue, Unit B, Brea, CA 92821

Alleva Medical®, extriCARE®, Vide®, Maxson® are Trademarks and/or Registered Trademarks manufactured by Alleva Medical Ltd., its subsidiaries, and/or affiliates in the United States, Europe, China and other countries. Copyright 2022.

All dimensions, weights, colours and contents may vary.